DiaLISS Control

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DiaMed GmbH
Pra Rond 23
1785 Cressier FR
Suisse

DiaLISS Control

English	B106701	03.13
for the control of the indirect Coombs test with DiaLISS		
Product-Identification: 12070		

INTRODUCTION

With the addition of a modified LISS-reagent, the incubation time for the indirect Coombstest can be reduced considerably while improving the sensitivity of the reaction (see instruction sheet 'DiaLISS'). To verify the reactivity of the LISS-reagent and the proper procedure of the test, the control should be included in each test-serie.

REAGENTS

'DiaLISS Control' is a human serum containing weak antibodies against the Rh antigenes D (RH1) and c (RH4). It reacts with all human red cells with the exception of such rare cases where both rhesus-antigenes D (RH1) and c (RH4) are missing. The result with 'DiaLISS Control' must be stronger than the result using isotonic saline solution. The actual result will vary according to the Rh phenotype of the red cells used.

Human serum contains Anti-c (RH4), Anti-D (RH1).
Preservative: < 0.1% NaN₃.

Caution: The source materials from which these products were manufactured, were found non reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.

ADDITIONAL REAGENTS REQUIRED

- 0.9% isotonic saline solution
- DiaLISS
- DiaClon Coombs-serum
- Coombs-control-cells

FURTHER MATERIALS REQUIRED

- Suspension tubes
- Tube rack
- Pipettes 50 µl, 1.0 ml
- Immunohaematological centrifuge
- Waterbath

PREPARATION OF BLOOD SAMPLE

Prepare a 3–5% red cell suspension in isotonic saline solution as follows:

1. Dispense 1.0 ml of isotonic saline solution into a clean glass tube.
2. Add 2 drops (100 µl) of whole blood or 1 drop (50 µl) of packed cells; mix gently.

TEST PROCEDURE

(see also instruction sheet 'DiaLISS')

1. Prepare 2 clean tubes.
2. Pipette 2 drops of the 'DiaLISS Control' serum into both tubes.
3. Add 1 drop of the 3–5% cell suspension* in isotonic saline solution into both tubes.
4. Add 4 drops of the 'DiaLISS' reagent into one of the tube.
5. Add 4 drops of isotonic saline solution into the other tube.
6. Mix and incubate both tubes for 5 minutes at 37 °C.
7. Wash contents of both tubes 3 times and decant supernatant.
8. Add 2 drops of 'DiaClon Coombs' serum to both tubes.
9. Mix well and centrifuge for 1 minute at 1000 rpm (RCF 125g) or for 20 seconds at 3400 rpm (RCF 1000g).
10. Gently resuspend the cell button and, over an indirect light source, examine macroscopically for agglutination.

* Remark: It is recommended to use a red cell suspension used in the test serie.

INTERPRETATION

The result with 'DiaLISS Control' must be stronger than the result using isotonic saline solution. The actual result will vary according to the Rh phenotype of the red cells used.

A negative reaction indicates an error in the test system; to verify, add 1 drop of Coombs-Control-cells to the negative tube, mix, centrifuge and observe agglutination.

- A positive reaction with Coombs-Control indicates an error of the LISS-test-system (too short incubation time, incorrect proportion of drops of serum – red cells – 'DiaLISS', LISS-reagent not reactive).
- A negative result with Coombs-Control indicates an error in the Coombstest system, with proteins (insufficient washing of the contents of a tube, contaminated wash solution, Coombs-serum was not added, Coombsserum is not reactive).

PRODUCTS

'DiaLISS Control'	1 x 10 ml	REF 106710
-------------------	---------------------	------------

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

Changes to the version 08.10 are shaded grey.

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